NMR OPERATION WITHOUT SAMPLE CHANGER

- 1. Wipe the exterior of the NMR sample tube using a clean Kimwipe to remove fingerprints and other deposits on the tube.
- 2. Carefully slide the sample into the *spectrometer* and proceed to sample analysis. The sample temperature will stabilize within 1-2 minutes, ensuring optimal performance while you are selecting/verifying the parameters for your NMR experiment:



Setting up basic parameters

• **Observe nucleus** – This button lets you choose the nucleus you wish to observe [1H, 13C, 19F, etc.]. Select the desired nucleus on the popup menu, followed by OK to confirm the selection.

1H	Nucleus	d.
13C		
		.1

• Lock nuclei - The lock nucleus is kept at the standard deuterium 2H for samples prepared in deuterated solvents. The OK confirms the selection.

2H	
1H	

• **Solvent** - Click the solvent button for a popup menu listing NMR solvents. Select the appropriate solvent, followed by OK to confirm the selection.

	Solvent	
D2O		
DMSO-d6		
Chloroform-d		
Methanol-d4		
Acetone-d6		
Acetonitrile-d3		
Benzene-d6		
TFA-d		
Ethanol-d6		
THF-d8		
Calibrate Solvent	Cancel	OK

• **Experiment** – Upon hitting this button, a popup window will appear with various 1-D and 2-D experiment options available. Select the appropriate experiment and the OK button to confirm.

Experi	menc
1D	
COSY	
JRES	
HSQC	
HMBC	
TOCSY	
Kinetics	
T1	
T2 - Hahn Ed	cho
T2 - CPMG	
Canad	04

• **Number of scans** - Upon hitting this button, a popup window with a list of the number of scans [in multiples of 4] can be specified depending on the desired spectral output. For a 1-H proton NMR 4 or 16 scans should give you a decent output.

Custom	
Custom:	
1	
1	
4	
16	
64	
256	
1024	
4096	
Cancel	Ok

NOTE: The instructor or technician will preset experiment settings/acquisition parameters depending on the experiment and the sample(s) being analyzed.

Once the basic parameters are set, hit the **GO** button to initiate the experiment.



After the first scan, the *Fourier-transformed* spectrum will be displayed and updated with signal-to-noise improvements with every additional 4 scans.



Saving the file

Upon hitting the disk icon on the screen, several saving options will appear. Click the export option to save the file to a USB key (thumb drive). Use a unique filename for each group. The default file format will be the *file*.dx format. Click the eject button to remove the thumb drive safely.



Remove your sample tube after the data has been saved.

On your lab computer, open the web browser, navigate to - <u>www.nmrium.org/nmrium</u>, and drag your xxxxxxx.dx file into the browser to begin spectral processing.

NMR SPECTRUM PROCESSING WITH NMRIUM

- 1. Open your web browser and navigate to <u>www.nmrium.org</u>
- 2. Drag your NMR file (xxxxxx.dx) into the browser to see the raw FID: (Skip to Step 7, if you already see a spectrum.)
- **A** 3. Click the **Apodization** button to reduce the noise. Set Line Broadening to 0.5 and click the Apply button.
- 4. Click the Zero Filling button to improve the resolution. Increase the Size to 32K and click the Apply button.
- K 5. Click the Fourier Transform button to convert the FID to a
- 6. Click the **Phase Correction** button. Change from Manual to Automatic and click Apply.
- J. Click **Baseline Correction** and click Apply to straighten the baseline.
- 3. Click the **Peak Picking** button and click Apply to display the peak positions.
 - 9. Click on the number above the TMS calibration signal and change it to "0".
- 10.Click the Integration button.
 - While holding the SHIFT key on the keyboard, click and drag the mouse pointer across the first signal. This will define the first integral region. Repeat for the other signals in your spectrum.



+ 11. Export your finished spectrum as a PNG image.





